


It's Neo-Hydro (Locking) /

Steps for use

These instructions are to be used as a reference guide only. Please refer to the full instructions for use prior to using this product.



Step 1.

Place the pigtail straightening device over the pigtail.

Step 2.

Feed the cannula (blue hub) through catheter and ensure the luer lock is secured.

Step 3.

Remove the pigtail straightening device.*

Step 4.

Remove the paper holder, advance the stylet and screw the green hub to the cannula.

Step 5.

Ensure the stylet tip is exposed beyond the catheter. The catheter is now ready to be inserted.

Step 6.

Unscrew the green luer lock hub and withdraw the stylet.

Step 7.

Attach a luer lock syringe to the clear plastic hub and aspirate fluid.**

Step 8.

Hold the catheter in place & unscrew the blue luer lock cannula hub. Keep the metal cannula still & slowly advance the catheter until you are satisfied with its position (arrow 1). Now withdraw the metal cannula (arrow 2).

Step 9.

Pull back the blue horse shoe clip and wrap the safety string around the inset within the white mould.

Step 10.

Place the horse shoe clip over the white mould. Blue plastic handle can be detached from the horse shoe clip & discarded.

* Please Note: To reduce risk of needle stick injury the pigtail straightening device can be removed by sliding it to the base of the catheter & peeling away. Dipping or wiping the Neo-Hydro in saline solution will activate the hydrophilic coating & facilitate smooth insertion.

** Please Note: If catheter placement is not correct at this stage the stylet can be re-introduced into the cannula and the system repositioned.