


## Kit contents:

- 28F delivery system with pre-loaded stent
- Balloon valve
- 50ml syringe
- 0.035 guidewire



## Ordering Information


| Product               | Product Code     | Nominal Diameter (mm) | Nominal Length |
|-----------------------|------------------|-----------------------|----------------|
| Danis Stent           | 019-08S-25-135-8 | 25                    | 135mm          |
| ELLA Extractor System | 015-01-28-1000   | 9.4                   | 100cm          |

## Additional information:

- Delivery system 28F / 20F, active length 60cm
- Retrieval of Danis requires use of the ELLA extractor system
- ELLA extractor assists the safe removal of Danis with minimal risk of trauma

## Related Products

| Product Code | Description               | Length | Tip           |
|--------------|---------------------------|--------|---------------|
| 3518004A     | Selectiva guidewire 0.35" | 180cm  | Angled, stiff |
| 353004A      | Selectiva guidewire 0.35" | 300cm  | Angled, stiff |



Manufactured by



References:

1. Wright GI, Lewis H, Hogan B, Burroughs, A, Patch D, O'Beirne J, 'A self-expanding metal stent for complicated variceal hemorrhage' (2010) Gastrointestinal Endoscopy
2. Escorsell A, et al., 'Esophageal Balloon Tamponade Versus Esophageal Stent in Controlling Acute Refractory Variceal Bleeding: A Multicenter Randomized Controlled Trial' (2016) Hepatology, (Vol 63, No.6)

Distributed by:

UK Medical  
Albreda House  
Lydgate Lane  
Sheffield, S10 5FH

T: +44 (0) 114 268 8880  
F: +44 (0) 114 268 8881  
sales@ukmedical.com  
www.ukmedical.com

