

UKMEDICAL®

It's Interventional.

*It's the Osteo-site®
vertebral bone
biopsy range!*



hello@ukmedical.com
ukmedical.com

It's about procedural choice.../

The Osteo-site® range of bone biopsy needles provides procedural flexibility and clinical choice.


Designed to ensure accurate placement and high quality sample capture, Osteo-site® devices are easy to use and offer a complete solution for your vertebral access procedures.

Product benefits:

- ***Controlled needle advancement***
- ***Precise visualisation for accurate needle placement***
- ***Quality sample capture***
- ***Available in a choice of tip designs to facilitate optimum placement***
- ***Range of gauges and lengths. Option for coaxial needle for multiple sampling capability***
- ***Designed in collaboration with Dr Kieran Murphy***
- ***Suitable for vertebral body biopsy and vertebroplasty procedures***
- ***Kit design minimises clinician radiation exposure during procedure***

Additional information?

If you would like further assistance with using our Vertebral Bone Biopsy Needles, or need advice on how to order, please call us on: **0114 268 8880**



Product features:



Ergonomic handle design provides tactile feedback

Radiolucent handle provides unobstructed view during introduction into bone



Option of 30° side-beveled stylet for minor redirection of needle during introduction

Option of sharpened diamond tip stylet provides optimum purchase

Centimetre depth markers


Murphy M1M™

Match ground stylet and cannula design

Coaxial option also available

Murphy M2™

Sharpened cutting cannula



Radiolucent needle holder




Diamond tip coaxial bone biopsy system

Murphy M2-S™

Ordering Information

Product Code	Size	Description	Global Product Number	Kit Contents
DBBN-11-10.0-M1M	11g x 10cm	Murphy M1M™ Side Bevel	G13761	
DBBN-11-15.0-M1M	11g x 15cm	Murphy M1M™ Side Bevel	G26949	
DBBN-13-10.0-M1M	13g x 10cm	Murphy M1M™ Side Bevel	G13762	
DBBN-13-15.0-M1M	13g x 15cm	Murphy M1M™ Side Bevel	G26950	
DBBN-11-10.0-M2	11g x 10cm	Murphy M2™-Diamond Tip	G13018	
DBBN-11-15.0-M2	11g x 15cm	Murphy M2™-Diamond Tip	G13020	
DBBN-13-10.0-M2	13g x 10cm	Murphy M2™-Diamond Tip	G13019	
DBBN-13-15.0-M2	13g x 15cm	Murphy™ M2-Diamond Tip	G13021	

Product Code	Size Outer Cannula	Size Inner Biopsy Needle	Description	Global Product Number	Kit Contents
DBBN-11-10.0-M2-S	11g x 10cm	14g x 17.1cm	Murphy M2-S™ Multi-Purpose Set	G36211	
DBBN-11-15.0-M2-S	11g x 15cm	14g x 22.1cm	Murphy M2-S™ Multi-Purpose Set	G36406	
DBBN-13-10.0-M2-S	13g x 10cm	16g x 17.1cm	Murphy M2-S™ Multi-Purpose Set	G36407	
DBBN-13-15.0-M2-S	13g x 15cm	16g x 22.1cm	Murphy M2-S™ Multi-Purpose Set	G36409	

Please see our vertebroplasty solutions brochure for more information.



Needle

+




Cement

+



Mixer

+



Injector

References:

* Please contact us for stock availability

** Please note that other needle types are available on request. These include: Cervical bone biopsy needle set; Diamond tip cervical bone biopsy needle set; Ratchet bone biopsy needle set; Coaxial curved bone biopsy needle set.
Please contact us for further information.

Manufactured by



Distributed by:

UK Medical
Albreda House
Lydgate Lane
Sheffield, S10 5FH

T: +44 (0) 114 268 8880
F: +44 (0) 114 268 8881
sales@ukmedical.com
www.ukmedical.com

