

UKMEDICAL®

It's Interventional.

*It's the Osteo-site®
vertebral bone
biopsy range!*



hello@ukmedical.com
ukmedical.com

It's about procedural choice.../

The Osteo-site® range of bone biopsy needles provides procedural flexibility and clinical choice.


Designed to ensure accurate placement and high quality sample capture, Osteo-site® devices are easy to use and offer a complete solution for your vertebral access procedures.

Product benefits:

- ***Controlled needle advancement***
- ***Precise visualisation for accurate needle placement***
- ***Quality sample capture***
- ***Available in a choice of tip designs to facilitate optimum placement***
- ***Range of gauges and lengths. Option for coaxial needle for multiple sampling capability***
- ***Designed in collaboration with Dr Kieran Murphy***
- ***Suitable for vertebral body biopsy and vertebroplasty procedures***
- ***Kit design minimises clinician radiation exposure during procedure***

Additional information?

If you would like further assistance with using our Vertebral Bone Biopsy Needles, or need advice on how to order, please call us on: **0114 268 8880**



Product features:



Ergonomic handle design provides tactile feedback

Radiolucent handle provides unobstructed view during introduction into bone



Option of 30° side-beveled stylet for minor redirection of needle during introduction

Option of sharpened diamond tip stylet provides optimum purchase


Centimetre depth markers

Match ground stylet and cannula design
Coaxial option also available

Murphy M1M™

Murphy M2™

Sharpened cutting cannula



Radiolucent needle holder



Diamond tip coaxial bone biopsy system

Murphy M2-S™



Ordering Information

Product Code	Size	Description	Global Product Number	Kit Contents
DBBN-11-10.0-M1M	11g x 10cm	Murphy M1M™ Side Bevel	G13761	
DBBN-11-15.0-M1M	11g x 15cm	Murphy M1M™ Side Bevel	G26949	
DBBN-13-10.0-M1M	13g x 10cm	Murphy M1M™ Side Bevel	G13762	
DBBN-13-15.0-M1M	13g x 15cm	Murphy M1M™ Side Bevel	G26950	
DBBN-11-10.0-M2	11g x 10cm	Murphy M2™-Diamond Tip	G13018	
DBBN-11-15.0-M2	11g x 15cm	Murphy M2™-Diamond Tip	G13020	
DBBN-13-10.0-M2	13g x 10cm	Murphy M2™-Diamond Tip	G13019	
DBBN-13-15.0-M2	13g x 15cm	Murphy™ M2-Diamond Tip	G13021	

Product Code	Size Outer Cannula	Size Inner Biopsy Needle	Description	Global Product Number	Kit Contents
DBBN-11-10.0-M2-S	11g x 10cm	14g x 17.1cm	Murphy M2-S™ Multi-Purpose Set	G36211	
DBBN-11-15.0-M2-S	11g x 15cm	14g x 22.1cm	Murphy M2-S™ Multi-Purpose Set	G36406	
DBBN-13-10.0-M2-S	13g x 10cm	16g x 17.1cm	Murphy M2-S™ Multi-Purpose Set	G36407	
DBBN-13-15.0-M2-S	13g x 15cm	16g x 22.1cm	Murphy M2-S™ Multi-Purpose Set	G36409	

Please see our vertebroplasty solutions brochure for more information.



Needle

+




Cement

+



Mixer

+



Injector

References:

* Please contact us for stock availability

** Please note that other needle types are available on request. These include: Cervical bone biopsy needle set; Diamond tip cervical bone biopsy needle set; Ratchet bone biopsy needle set; Coaxial curved bone biopsy needle set.

Please contact us for further information.

Manufactured by



Distributed by:

UK Medical
Albreda House
Lydgate Lane
Sheffield, S10 5FH

T: +44 (0) 114 268 8880
F: +44 (0) 114 268 8881
sales@ukmedical.com
www.ukmedical.com

