

UKMEDICAL®

It's Interventional.

*It's
BiopSafe® /*



hello@ukmedical.com
ukmedical.com

It's about removing a recognised risk... /

One of the most common problems with biopsy handling is the risk of being exposed to formalin and formaldehyde either through touch or inhalation.

Recent EU regulations reclassified formaldehyde as a category 1b carcinogen and category 2 mutagen.

Suited to a wide range of biopsy procedures, BiopSafe offers the same quality sample fixation you're used to, but with additional safety.

Product benefits:

- ***Safe biopsy handling with no exposure to formalin or formaldehyde via touch or inhalation***
- ***Gold standard sample preparation***
- ***Reliable and trusted***
- ***Ready in seconds***
- ***No sensory irritation***
- ***Optimal formalin volumes minimises waste and expensive disposal costs***

DANGER

**FORMALDEHYDE
IRRITANT & POTENTIAL
CANCER HAZARD**

Additional information?

If you would like further assistance with using BiopSafe, or need advice on how to order, please call us on: **0114 268 8880**

Did you know...

As soon as you open the lid on a prefilled formalin pot you risk being exposed to formaldehyde?

"Formaldehyde is carcinogenic to humans" - International Agency for Research on Cancer¹



"A known human carcinogen" - National Toxicology Program²

"May cause leukaemia, particularly myeloid leukaemia in humans" - National Cancer Institute³

Product features:

The secret is in the lid!

The BiopSafe is a small, easy to use sample pot with a formalin capsule built in to the lid



The formalin is released by simply applying pressure to the button on the lid

Everything is done in the sample pot



Optimised formalin volumes minimise waste and expensive disposal costs

10% neutral buffered formalin

The lid and pot form a hermetic seal, ensuring there's no risk of leaking or exposure



Anti cross-threading design ensures hermetic seal



Available in 20ml and 60ml sizes to cover a wide range of tissue types and volumes

Ordering Information

Product Code	Description	Volume (ml)	Suited for sample sizes up to
XL-3178-200-00	BiopSafe Biopsy container 20ml	20	1 cm ³
XL-3178-600-00	BiopSafe Biopsy container 60ml	60	3 cm ³



Suitable for use with biopsy cassettes.

Related Items

Our range of Achieve, Temno and TruCut soft tissue biopsy needles are available on request alongside our Jamshidi bone marrow biopsy needles.

Support

- Training packs including IFU's, presentations, videos and scientific papers available upon request
- Individual and departmental training sessions available
- Get in touch via: 0114 268 8880 or hello@ukmedical.com

Manufactured by



References:

1. International Agency for Research on Cancer (2009). IARC monographs on the evaluation of carcinogenic risks to humans v. 100F. pg 430.
2. NTP (National Toxicology Program). 2016. Report on Carcinogens, Fourteenth Edition.; Research Triangle Park, NC: U.S. Department of Health and Human Services, Public Health Service. <http://ntp.niehs.nih.gov/go/roc14/> (HYPERLINK "/pubhealth/roc/14th_endnote.xml"EndNote XML).
3. National Cancer Institute. Formaldehyde and Cancer Risk. (June 2011) <https://www.cancer.gov/about-cancer/causes-prevention/risk/substances/formaldehyde/formaldehyde-fact-sheet#r2> Last accessed 18/08/2017.

Distributed by:

UK Medical
Albreda House
Lydgate Lane
Sheffield, S10 5FH

T: +44 (0) 114 268 8880
F: +44 (0) 114 268 8881
sales@ukmedical.com
www.ukmedical.com

@UK_Medical

 **UKMEDICAL[®]**
It's Interventional.