

It's Epi-Guard/


Fixation Device for Epidural Catheters





Epi-Guard is a fixation device for epidural catheters. The device consists of a skin-side surface with a skin friendly adhesive combined with a special design.

Product benefits:

- *Reduces risk of accidental removal of catheter*
- *Safe anchoring and protection of catheter*
- *Allows inspection of puncture site and catheter position*
- *Can prevent kinking of catheter*
- *Can prevent contamination of puncture site from environment*
- *Easy replacement of anchoring device*



Product features



Safe anchoring of catheters

Epi-Guard is an anchoring device for epidural catheters. The device consists of a skin-side surface with a skin friendly adhesive combined with a specifically designed foam block that safely anchors epidural catheters.

The hydrocolloid dressing will absorb body fluids or any anesthetic given prior to insertion of the epidural catheter which can flow back from the needle puncture site.

The soft foam support block can prevent kinking of the catheter, even if the patient should lie on it.

For the sake of the patient

Epidural catheters are typically secured inadequately by PU-film or other adhesive dressings and accidental displacement of the catheter is frequently seen. Not only will this cause the patient to experience a period without sufficient anesthesia, but a very uncomfortable insertion of a new catheter is also required.

When Epi-Guard is being used, accidental displacement of the epidural catheter will be rare. The hydrocolloid carrier has been carefully tested for skin irritation and the risk of skin reactions is minimal.

The Epi-Guard gives the patient great freedom to move around without worrying about the catheter.

Easy to use

The Epi-Guard is easy and fast to apply to the patient's skin and results in a standardised and safe anchoring of the epidural catheter.

Ordering Information

Description	Design type	Product code	Box quantity
Epi-Guard fixation device	Oval - Anti-Kink - Hydrocolloid	AD-8170	25

Distributed by:

*UK Medical
Albreda House
Lydgate Lane
Sheffield, S10 5FH*

*T: +44 (0) 114 268 8880
F: +44 (0) 114 268 8881
sales@ukmedical.com
www.ukmedical.com*

 *@UK_Medical*

 **UKMEDICAL[®]**
It's Interventional.