

UKMEDICAL®

It's Interventional.

***It's the UK Medical
Nephrostomy Service***

District Nurse Information



Product information



How do I place an order?

- Your patient will need replacement drainage leg bags, these are recommended to be changed weekly. Please see NHS SC code below for order details. Note, the 800ml leg bag should be emptied in accordance with your normal regime, but does not require an overnight bag.
- The dressing (catheter fixation device) will need to be changed on a weekly basis. Please select preferred device below and right.

Product Code	Description	Contents	NHS Supply Chain Code
DB800C-MS-V	Leg bag 800ml	800ml leg bag with male luer lock and velcro strap	FRZ3891
AD-8150	Drain Guard	Catheter fixation device	FRJ8133
BT-CS-0512	Drain-Lok	Catheter fixation device	FRZ3892

How do I place an order direct?

- To ensure continuity of care and speed of delivery, our nephrostomy equipment is not available on prescription. These can be ordered directly from UK Medical.
- Direct orders with UK Medical can be placed by fax, telephone and email:
 - Fax: 0114 268 8881
 - Telephone: 0114 268 8880
 - Email: sales@ukmedical.com
- Orders placed before 3:30pm will be despatched the same day, for next day delivery directly to your preferred address.
- If you have an urgent order after 3:30pm, please call our Nephrostomy Support Team on 0114 268 8880.
- Our nephrostomy equipment can also be ordered through NHS Supply Chain using the codes above (right).
- UK Medical cannot guarantee delivery times on orders placed via NHS Supply Chain.

Additional information?

- If you would like further assistance with using our nephrostomy equipment, or need advice on how to order, please contact please contact your UK Medical Area Business Manager
- Alternatively contact our Nephrostomy Support Team on **0114 268 8880** or email customersupport@ukmedical.com

Catheter fixation

Drain Guard: AD-8150



Perforated cover for easy placement of anchoring device.

Transparent hydrocolloid cover allows monitoring of puncture site.

Smooth non-woven cover with strong adhesion for secure catheter fixation.

PE-Foam support.

Thin, comfortable hydrocolloid carrier.



Drain-Lok: BT-CS-0512

In PET/Rayon.


TPU Cushion.

Special design silicone support.

Hydrocolloid adhesive with an outer impermeable film.

Liner removal tab.

Extra-strong adhesive ensures fixation on the TPU Cushion.



For more information about Neo Hydro or other Bioteq products, please contact our **Customer Support Team on 0114 268 8880, email: customersupport@ukmedical.com or visit: www.ukmedical.com**

Manufactured by

BIOTEQ[®]

Distributed by

*UK Medical
Albreda House
Lydgate Lane
Sheffield, S10 5FH*

*+44 (0) 114 268 8880
hello@ukmedical.com
www.ukmedical.com*

2021-02-12-Nephrostomy Service District A5 Blank



UKMEDICAL[®]

It's Interventional.